# 1 FILTRATION

Ashless Filter papers for qualitative analysis

General Purpose Filter papers for qualitative analysis

Ashless Filter papers for quantitative analysis

Ashless Hardened Filter papers for quantitative analysis

Hardened Low Ash Filter papers for quantitative analysis



## **CHM®** Ashless **Filter papers for qualitative analysis**

These filter papers are use for qualitative analysis. Qualitative filters are made of refined pulp and linters with an alpha-cellulose content of virtually 100%, which gives them a number of diverse filtration properties. The ash content of this filter paper is not reduced by post-treatment. Qualitative filter papers are available as sheets, discs and folded filters.

### F1001 GRADE - Medium Filtration

The most widely used filter paper in the CHM range. Medium retention and flow rate. This grade covers a wide range of laboratory applications and is frequently used for clarifying liquids. Traditionally the grade is used in qualitative analytical separations for routine laboratory work, rapid filtration of fine precipitates such as lead sulphate, calcium oxalate (hot) and calcium carbonate. In agriculture, it is used for soil analysis and seed testing procedures. In the food industry, Grade F1001 is used for numerous routine techniques to separate solid foodstuffs from associated liquid or extracting liquid and is widely used in education for teaching simple qualitative analytical separations. In air pollution monitoring, using circles or rolls, atmospheric dust is collected from an airflow and the stain-intensity measured photometrically. For gas detection, the paper is impregnated with a chromogenic reagent and colour formation quantified by optical reflectance. F1001 circle sizes available from 25 to 650mm diameter. Sheets and strips also available.

### F1002 GRADE - Medium Slow Filtration

Slightly more retentive than No. F1001 and therefore with a moderate to slow filtration speed. More retentive and absorbent than No. F1001. In addition to general filtration this grade F1002 is used, for example, monitoring specific contaminants in the atmosphere, filtration of fine precipitates, soil testing, often used as folded filter in an analytical funnel.



Technical Specifications							
GRADE	Properties	Weight g/m²	Thickness µm	Retention range µm	Ash content %		
F1001 F1002 F1003 F1004 F1005 F1006	Medium Medium-Slow Medium-Slow/Thick Very fast Very Slow Slow	85 85 200 85 85 80	190 180 320 210 170 150	10 - 13 7 - 9 5 - 7 20-25 3 - 5 2 - 4	<0.06 <0.06 <0.06 <0.06 <0.06 <0.06		

GRADE Dia. (mm)	F1001	F1002 Standard Gra	F1003 ades Circles (10	F1004 00/box)	F1005	F1006
42,5 47 55 70 90 110 125 150 185 240 320	F1001-042 F1001-047 F1001-055 F1001-070 F1001-090 F1001-110 F1001-125 F1001-150 F1001-185 F1001-240 F1001-320	F1002-042 F1002-047 F1002-055 F1002-070 F1002-090 F1002-110 F1002-125 F1002-150 F1002-185 F1002-240 F1002-320	F1003-042 F1003-047 F1003-055 F1003-070 F1003-090 F1003-110 F1003-125 F1003-150 F1003-185 F1003-240 F1003-320	F1004-042 F1004-047 F1004-055 F1004-070 F1004-090 F1004-110 F1004-125 F1004-150 F1004-185 F1004-240 F1004-320	F1005-042 F1005-047 F1005-055 F1005-070 F1005-090 F1005-110 F1005-125 F1005-150 F1005-185 F1005-240 F1005-320	F1006-042 F1006-047 F1006-055 F1006-070 F1006-110 F1006-125 F1006-150 F1006-185 F1006-240 F1006-320
GRADE Size (mm)	F1001	F1002 Standard Gra	F1003 ades Sheets (10	F1004 00/pack)	F1005	F1006
460x570	F1001-460570	F1002-460570	F1003-460570	F1004-460570	F1005-460570	F1006-460

### F1003 GRADE - Medium Slow Filtration (Thick)

Medium to low rate of filtration with double the thickness comparing with CHM Grade F1001. Fine particle retention and excellent loading capacity. The extra thickness gives increased wet strength and allows a higher solute loading. Preferably used for liquids hard to clarify, essences, oils, tinctures, particularly useful for use in Buechner-funnels.

### F1004 GRADE - Fast Filtration

Very high rate of filtration with excellent retention of coarse precipitates such as metal hydroxides and sulphides, or gelatinous substances. Preferably used for as a rapid filter for various organic metal precipitates, routine cleanup of biological fluids, food industry analysis, air pollution monitoring (high rates and the fine particles collection is not critical).

### F1005 GRADE - Slow filtration

Lowest rate of filtration in the CHM qualitative range and maximum degree of fine particle filtration. Most efficient filtration of smallest particles. Preferably used as clarifying filter for cloudy suspensions and for water and soil analysis. Particularly in diffucult filtration conditions and extra fine-grained precipitates, barium sulphate, cuprous oxide, often specified for clarification of wine.

### F1006 GRADE - Slow Filtration.

Similar particle retention as Grade F1005 with higher filtration speed. Often used for boiler water analysis.



## **CHM®** General Purpose **Filter papers for qualitative analysis**

These general purpose filters have a high Wet Strengthened. They are made of high-purity cotton linters and other vegetable fibres. These filter papers have fast or very fast filtration rates, and are particularly useful in filtering coarse precipitates or relatively straightforward substances. These filters are not recommended for Kjeldahl.

### F1091 GRADE - Very fast filtration

Creped surface. Thick Creped filter paper with medium flow rate. For general laboratory use in less-critical analyses. Used around the world in laboratories to assay sugar cane or beet. The fruit is mashed and further analyzed according to the aluminium sulfur method.

### F1093 GRADE - Fast filtration

Smooth Grade F1093 is a General Purpose filter paper for qualitative analysis. This wet strengthened paper is used for general filtration and sample preparation in the food sector, sugar processing plants, hospitals, educational and research centres, colleges, universities and labs (with a very high usage and less critical analysis), etc.

### F1113 GRADE - ExtraFast filtration. Thick

High particle retention and extremely high loading capacity. Preferably used for filtration of gelatine, resins solutions and other viscous liquids, such as syrups, oils, essences and fats. The folded format enables greater volumes to be dealt with at atmospheric pressures.







### **Technical Specifications** GRADE Properties Weight Thickness Retention range Ash content F1091 Very Fast. Crepped. Thick 34-42 64 165 n/a F1093 80 180 35-40 Very Fast. n/a F1113 Extra-Fast. Thick 470 60-65 n/a

GRADE	F1091	F1093	F1113
Dia. (mm)	Stai	ndard Grades Circles (100/box)	
42,5	F1091-042		
47	F1091-047		
50	F1091-050		
55	F1091-055	F1093-055	F1113-055
70	F1091-070	F1093-070	F1113-070
90	F1091-090	F1093-090	F1113-090
110	F1091-110	F1093-110	F1113-110
125	F1091-125	F1093-125	F1113-125
150	F1091-150	F1093-150	F1113-150
185	F1091-185	F1093-185	F1113-185
200	F1091-200	F1093-200	F1113-200
240	F1091-240	F1093-240	F1113-240
320	F1091-320	F1093-320	F1113-320
400	F1091-400	F1093-400	F1113-400
500	F1091-500	F1093-500	F1113-500
650	F1091-650	F1093-650	F1113-650
GRADE	F1091	F1093	F1113
Size (mm)	Sta	ndard Grades Sheets (100/box)	
460 x 570	F1091-460570	F1093-460570	F1113-460570



## **CHM®** Ashless **Filter papers for quantitative analysis**

These CHM filter papers are used for quantitative analyses and designed for preparation of samples and gravimetric analysis. They are made of refined pulp and linters with virtually 100% of alpha-cellulose content. These filter papers are guaranteed free of possible residual acids used in some production methods. Another main feature is the extremely low percentage of ash content (maximum ash content of 0.007%). Chemiton offers you two groups of ashless filter papers for quantitative analyses. They are suitable for Buechner funnels and for filtration under pressure.

### F2040 GRADE - Medium Slow filtration

The classic general purpose ashless filter paper with a medium-to-slow filtering rate. Suitable for typical applications include gravimetric analysis for numerous components and for pre-filtrations of all kinds. Preferably used for as a primary filter for separating solid matter from aqueous extracts, various tests for fat and oil in water, in general soil analysis, quantitative determination of sediments in milk, as an analytical grade clean-up filter for solutions prior to AA spectro- photometry. Suitable for finer precipitates, such as hot barium sulphate.

### F2041 GRADE - Fast filtration

Fast ashless filter paper in the CHM quantitative range. It is particularly suitable for analytical procedures and tests involving large particles or gelatinous precipitates (e.g metal hydroxides and sulphides). Also used in metal (Pb) tests in water testing analysis, quantitative air pollution analysis, food industry, paper industry, etc.

Order Information						
GRADE Dia. (mm)	F2040	0	F2042 Grades Circles	F2043 s (100/box)	F2044	F2045
12,5 47 55 70 90 110 125 150 185 240	F2040-047 F2040-055 F2040-070 F2040-090 F2040-110 F2040-125 F2040-150 F2040-185 F2040-240	F2041-047 F2041-055 F2041-070 F2041-090 F2041-110 F2041-125 F2041-150 F2041-185 F2041-240	F2042-047 F2042-055 F2042-070 F2042-090 F2042-110 F2042-125 F2042-150 F2042-185 F2042-240	F2043-012 (*) F2043-047 F2043-055 F2043-070 F2043-110 F2043-125 F2043-150 F2043-185 F2043-240	F2044-047 F2044-055 F2044-070 F2044-090 F2044-110 F2044-125 F2044-150 F2044-185 F2044-240	F2045-047 F2045-055 F2045-070 F2045-090 F2045-110 F2045-125 F2045-150 F2045-185 F2045-240
Size (mm) Ashless Grades Sheets (100/box)						
460 x 570	F2040-460570	F2041-460570	F2042-460570	F2043-460570	F2044-460570	F2045-460570
* Packs of 10	000 circles					

Techr	Technical Specifications							
GRADE	Properties	Weight g/m²	Thickness µm	Retention range µm	Ash content %			
F2040 F2041 F2042 F2043 F2044 F2045	Medium-Slow Fast Very Slow Medium Slow Very fast	80 80 100 80 80 80	170 190 160 180 160 210	7 - 9 20 - 25 2 - 3 14 - 17 2 - 4 25 - 30	<0.007 <0.007 <0.007 <0.007 <0.007 <0.01			

## F2042 GRADE – Very Slow filtration

An ashless world standard filter for critical gravimetric analysis. With a slow filtering rate and fine particle retention. Typical analytical precipitates such as cold barium sulphate, lead sulphate, zinc and nickel sulphides, etc. Using the methods of gravimetry, photometry, colorimetry and precipitation in drinking and waste water treatments it is suitable for testing for chemical elements and for some radioactive substances.

## F2043 GRADE - Medium Fast filtration

Ashless filter with medium filtration speed and good retention (between GRADE F2040 and GRADE F2041) of medium and thick particles. Suitable for Gravimetric measurements of gypsum/lime suspensions in power plants. F2043 GRADE is particularly applied in metallurgical industry laboratories for metal tests. Typical applications include foodstuffs analysis; soil analysis; particle collection in air pollution monitoring, COD and TOC determination, inorganic analysis in the construction, mining and steel industries, for the Blaine test in the cement industry for ascertaining how fine the cement is (standards UNE 80-112-91 and EN 196-6), and for carrying out other chemical analyses on cement.

### F2044 GRADE - Slow filtration

Thin version of No. F2042 but with higher flow rate (twice as fast as No. F2042). Very fine particles but with lower ash weight per sample.

### F2045 GRADE - Fast filtration

Filter paper of very high rate of filtration, wide-pored, soft, spongy structure, extremely low-ash content. For applications like: Food industry: determination of ash contents in foodstuffs, PCB determination in foodstuffs; Beverage industry: Processing (ashing) fruit juice samples for photometric determinations (e.g. phosphate); Environmental analysis: Determination of filterable substances and the residue on ignition (dry weight) according to the German Standard Methods for the examination of water, waste water and sludge (DIN 38 409 part 2)





## **CHM®** Ashless Hardened **Filter** papers for quantitative analysis

Hardened Ashless Filter papers are acid hardened which reduces the ash content to an extremely low level. This filters are produced by a complex, elaborate washing process under stringently controlled conditions. Firstly, acid washing is arranged. Then comes a series of washes in demineralised water, which increase the strength of the paper, thus making them particularly suitable for Buechner filter funnels and for a wide range of critical analytical filtration operations. Through this process, a maximum ash content of 0.007% is attained, which means that no contaminants are introduced when filtering, and also that full compliance with international standards on this subject is achieved.

Technical Specifications							
GRADE	Properties	Weight g/m²	Thickness µm	Retention range µm	Ash content %		
F2140 F2141 F2142	Medium Very fast Slow	80 80 80	160 170 150	7 - 12 20 - 25 2 - 4	0,006 0,006 0,006		



### F2140 GRADE - Medium Fast Filtration

Hardened ashless filter paper with medium retention and flow rate. Extremely strong and pure. With a hard surface is recommended for filtering medium-sized precipitates such as most metal sulphides. High chemical resistance. Used in the gravimetric analysis of metals in acid and slightly alkalinized solutions, pressure filtration

### F2141 GRADE - Fast Filtration

Hardened ashless filter paper with fast flow rate. Preferably used in filtration of coarse-flocculent and bulky precipitates (as aluminium, chromium or hydroxides of iron, bismuth, cobalt, sulphides of copper, various organic metal precipitates, etc.) and gelatinous precipitates in acid/alkali solutions during gravimetric analysis.

### F2142 GRADE - Slow Filtration

Hardened ashless filter paper with high retention and slow flow rate. High chemical resistance. Often used for filtering very fine precipitates and in gravimetric metal determinations.

Order Information							
GRADE Dia. (mm)	F2140 Hardei	F2141 ned Ashless (10	F2142 10/box)	F2050 Hardened	F2052 d Low Ash (100)	F2054 /box)	
42,5 47 55 70 90 110 125 150 185 200 240 320	F2140-070 F2140-090 F2140-110 F2140-125 F2140-150 F2140-185 	F2141-070 F2141-090 F2141-110 F2141-125 F2141-150 F2141-185 	F2142-070 F2142-090 F2142-110 F2142-125 F2142-150 F2142-185  F2142-240	F2050-042 F2050-047 F2050-055 F2050-070 F2050-090 F2050-110 F2050-125 F2050-150 F2050-185 F2050-200 F2050-240 F2050-241	F2052-042 F2052-047 F2052-055 F2052-070 F2052-090 F2052-110 F2052-125 F2052-150 F2052-185 F2052-200 F2052-240 F2052-241	F2054-042 F2054-047 F2054-055 F2054-070 F2054-090 F2054-110 F2054-125 F2054-150 F2054-185 F2054-200 F2054-240 F2054-241	
Size (mm)	Hardened Ashless (100/box) Hardened Low Ash (100/box)					/box)	
460 x 570				F2050-460570	F2052-460570	F2054-460570	



## **CHM®** Hardened Low Ash Filter papers for quantitative analysis

These filters, made from cotton linters fibre, are put through a washing process and treated with strong acids, and then baths in demineralized water, to produce high wet strength (makes them appropriate for filtering in low pressure or vacuum conditions) and chemical resistance (makes them suitable to work with acids or alkaline solutions in moderate concentrations). A very low ash-content filter with a 0.06% (the maximum ash contents of these filters is intermediate between CHM qualitative grades and ashless quantitative grades), very smooth surface, it makes easy to recover the whole of the precipitate after the filtration which is particularly indicated for Buchner filtrations.

### F2050 GRADE - Slow Filtration

It is the thinnest CHM filter with slow filtering rate, with excellent retention of very fine particles, such as barium sulphate, zinc sulphide, etc. Hardened and glazed surface makes this paper suitable for use in the electronical industry in carriers of electronical components or boards.

### F2052 GRADE - Medium Fast Filtration

General purpose hardened filter paper with medium-rate filtering, with good retention of medium particles, such as calcium oxalate and metal sulphides. Suitable for various tests on the intake of atmospheric pollution (sulphur oxides, ammonia gases, etc.) and also for the microbiological analysis of water.

They are used in fat extraction equipment as well in the oilseed and food industries, and in a large number of routine analytic procedures. Ask for our 3 pieces filter funnels.

### F2054 GRADE - Fast Filtration

The fastest filter paper in the range. Suitable for filtering coarse, gelatinous or dense liquids. Good load capacity.



12

Technical Specifications							
GRADE	Properties	Weight g/m²	Thickness µm	Retention range µm	Ash content %		
F2050 F2052 F2054	Slow Medium-Fast Fast	90 90 90	180 190 200	3 - 4 7 - 8 20 - 25	<0.06 <0.06 <0.06		



Glass Microfibre Filters

Quartz Microfibre Filters

pH indicator papers

Filter Paper in Reams

 $m{arphi} m{arphi} m{arphi} m{arphi}$  13